

**STANDARD OPERATION PROCEDURE**  
**Faculty of Biosciences, NMBU**

---

**Method name: Urea**

BIOVIT-nr.: Arb\_1012\_feed and muscle

---

**1. Introduction**

The method is intended for the determination of urea in feed and muscle.

The fish fed on diets that differed in levels of urea supplementation impact the incidence of salmon with ulcer, and subsequently its mortality that seems to relate to plasma osmolality. According to Rürvik et al. (2000) the dietary urea supplementation may lead to reduced development of skin ulcers and will probably increase the proportion of market size salmon of superior quality.

**Main instrument:** RX Daytona + (Randox Laboratories Ltd, UK). 55 Diamond Road, Crumlin, County Antrim, BT29 4QY, United Kingdom. Kit: Urea (UR8334).

Urea concentration can be measured photometrically on a RX Daytona + spectrophotometer. Reportable range: 0.50-62.0 mmol/L

**2. Reagents:**

Randox UREA Assay (R1 and R2)  
Clinical Chemistry Calibrator Level 2 CAL2350  
Clinical Chemistry Calibrator Level 3 CAL2351  
Quality control (Assayed Chemistry Premium Plus Level 3, Cat No-HE1532)

**3. Risk assessments:**

None.

**4. Recommended equipment:**

50 mL test tube (Greiner) with screw cap, and V bottom  
Erlenmeyer flask  
Filter paper (Whatman™)  
Double Spatula  
Tweezers  
Pipettes and pipette tips 1-5 mL  
Dispenser tips Combitips advanced® Eppendorf Quality, 50mL  
Bench centrifuge with a swing-out rotor  
Weighing scale

BIOVIT/NMBU						ARB
Prepared by Milena Bjelanovic	Approved by Hanne Kolsrud Hustoft	Valid from 092022	Revision	Replaced	Document name Arb1012 Urea_feed and muscle	Page 1-3

Vortex mixer  
RX Daytona + instrument

A control sample containing 5.1-6.0 % urea, should be run with each batch.  
**For feed analysis**, 1g ± 5 mg sample is needed. Degree of grinding: 0.5 mm.  
**For urea muscle analysis**, a 15-20 g ± 50 mg of homogenized sample is needed.

Feed samples must be stored from 2 to 8 °C.  
Muscle sample must be stored at -20 °C.

## 5. Job description

- a) If feed sample, weigh in about 1 g into 50 mL plastic tube w/ screw cap, or weigh in about 20 g if muscle sample.
- b) Make a note of the exact sample weight.
- c) Add Milli-Q water up to 50ml mark in the sample tube to the sample and shake it well for 10 sec on Vortex. Note the exact amount of added water to each sample.
- d) Centrifuge the samples for 15 min at 2500g and 10°C.
- e) Filter the sample and collect the supernatant in Erlenmeyer flask.
- f) Pipette 300 mL of the supernatant into Eppendorf tubes without caps.
- g) The samples are analyzed on a RX Daytona + fully automatic analyzer (See in MpsArb\_2020 Kliniske analyser Daytona +).

## 6. Calculation

% Urea= ((abs sample\*MW) \*0,05) / mg sample)/10  
abs. sample = urea absorbance read on spectrophotometer (mmol / l)  
60,06 = MW-molecular weight of urea (mg / mmol)  
0.05 = Amount of ml water added  
mg sample = weighed sample  
10 = the calculation in %

## 7. Literature:

<https://journals.physiology.org/doi/full/10.1152/advan.00027.2002>

<https://acutecaretesting.org/en/articles/urea-and-the-clinical-value-of-measuring-blood-urea-concentration>

[https://www.researchgate.net/publication/230107985\\_Urea\\_in\\_feeds\\_for\\_sea\\_water\\_farmed\\_Atlantic\\_salmon\\_Effect\\_on\\_growth\\_carcass\\_quality\\_and\\_outbreaks\\_of\\_winter\\_ulcer](https://www.researchgate.net/publication/230107985_Urea_in_feeds_for_sea_water_farmed_Atlantic_salmon_Effect_on_growth_carcass_quality_and_outbreaks_of_winter_ulcer)

BIOVIT/NMBU						ARB
Prepared by Milena Bjelanovic	Approved by Hanne Kolsrud Hustoft	Valid from 092022	Revision	Replaced	Document name Arb1012 Urea_feed and muscle	Page 2-3

BIOVIT/NMBU						ARB
Prepared by Milena Bjelanovic	Approved by Hanne Kolsrud Hustoft	Valid from 092022	Revision	Replaced	Document name Arb1012 Urea_feed and muscle	Page 3-3